

Comparison of AMS 700™ CX and Coloplast Titan® Inflatable Penile Prosthesis Cylinders Subjected to In-Vitro Cyclic Buckling

Charles Pritchard, MD
Geoff Daniels, BSME

Poster Presentation #111
Presented by Dr. Charles Pritchard
Sexual Medicine Society Meeting
October 2008
Toronto, CAN

ABSTRACT

Background: Parylene coating has been proven to increase the durability of AMS 700 series penile prosthesis cylinders. However, aneurysmal failures in parylene-coated AMS 700 cylinders during “woman on top” intercourse have been reported. It is surmised that buckling of the cylinder results in this particular mode of failure. We attempt to compare the resistance of AMS 700 and Coloplast Titan cylinders to this type of failure.

Methods: Five (1 non-parylene-coated and 4 parylene-coated) 18cm AMS 700 CX and 2 18cm Coloplast Titan cylinders were compared. Each cylinder was fixed proximally and distally, with the distal tip placed in a 10 degree angled fixture to simulate the vaginal canal angle. The cylinders were inflated to 15 PSI and placed in a 98.6°F saline bath. Axial loading was created using an Instron electrodynamic test system, resulting in buckling of the cylinder. Cylinders were examined every 100,000 cycles for aneurysms and degradation.

Results: In the non-parylene coated AMS 700 cylinder, no aneurysm formed, but the outer silicone layer failed at 600,000 cycles. The 4 parylene-coated cylinders developed aneurysmal defects between 1.1 and 2.05 million cycles, 2 also with catastrophic failures. The aneurysms occurred in areas of separation of the fibers in the fabric layer. The 2 Titan cylinders have been cycled 3.5 and 6.5 million times without signs of compromise in the integrity of the cylinder wall.

Conclusion: The durability of AMS 700 cylinders has been improved by the addition of parylene, but they may be prone to aneurysm formation due to buckling of the cylinders. Initial in-vitro studies suggest that Coloplast cylinders are inherently more resistant to this type of failure. Further testing will be required to fully define these differences and their clinical significance.

COLOPLAST KEY TAKEAWAYS

- In the non-parylene coated AMS 700 cylinder, no aneurysm formed, but the outer silicone layer failed at 600,000 cycles.
- The 4 parylene-coated cylinders developed aneurysmal defects between 1.1 and 2.05 million cycles, 2 also with catastrophic failures. The aneurysms occurred in areas of separation of the fibers in the fabric layer.
- The 2 Titan cylinders have been cycled 3.5 and 6.5 million times without signs of compromise in the integrity of the cylinder wall.
- Initial in-vitro studies suggest that Coloplast cylinders are inherently more resistant to this type of aneurysmal failure.

BRIEF STATEMENT

Indications: The Titan® and Titan Touch Penile Prosthesis is indicated for male patients suffering from erectile dysfunction (impotence) who are considered to be candidates for implantation of a penile prosthesis.

Contraindications: The Titan, and Titan Touch Inflatable Penile Prosthesis is contraindicated in patients who have one or more of the following: Patients with an active infection present anywhere in the body, especially urinary tract or genital infection. Patients with a documented sensitivity to silicone. Patients with unresolved problems affecting urination, such as an elevated residual urine volume secondary to bladder outlet obstruction or neurogenic bladder. Patients unwilling to undergo any further surgery for device revision.

Warnings: Implantation of the device may make latent natural erections, as well as other interventional treatment options, impossible. Men with diabetes or spinal cord injuries, as well as immunocompromised patients, may have an increased risk of infection associated with a prosthesis. Failure to evaluate and promptly treat erosion may result in a substantial worsening of the condition, leading to infection and loss of tissue. Implantation of a penile prosthesis may result in penile shortening, curvature or scarring. Pre-existing abdominal or penile scarring or contracture may make surgical implantation more complicated or impractical.

Precautions: Surgeons implanting penile prostheses should be familiar with the currently available techniques for measuring the patient, determining implant size, and performing the surgery. Removal of an implanted prosthesis without timely reimplantation of a new prosthesis may complicate subsequent reimplantation or make it impossible. A thorough preoperative consultation should include a discussion between the patient and physician of all available treatment options and their risks and benefits. The prosthesis should not be implanted in patients who lack the manual dexterity or strength necessary to operate the device. The device may be used in the presence of Peyronie's Disease. If the manual modeling technique is to be utilized, see the Surgical Protocol for more information.

Potential Complications: Scrotal swelling, auto-inflation, discomfort, angulation/curvature, edema, device malfunction, chronic pain, difficulty with ejaculation, transient urinary retention, fever, migration, patient dissatisfaction, site infection, deflation, hematoma/seroma, wound leakage, bleeding, delayed wound healing, phimosis, sensory loss, cylinder aneurysm, fibrous capsule formation, over/under inflation, erosion scrotal erythema, genital change, wound infection, and inguinal hernia.

The information provided is not comprehensive with regard to product risks. For a comprehensive listing of indications, contraindications, warnings and precautions refer to the product's Instructions for Use. Alternatively, you may contact a Coloplast representative at 1-800-258-3476 and/or visit the company website at www.coloplast.com.

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

[Ostomy Care / Continence Care / Wound & Skin Care / Interventional Urology](#)

Coloplast Corp. Minneapolis, MN 55411 / Interventional Urology Surgical Support 1-800-258-3476
www.coloplast.us The Coloplast logo is a registered trademark of Coloplast A/S.

© 2022 Coloplast Corp. All rights reserved.



PM-06072 02.22