


Genesis[®]

Malleable Penile Prosthesis

The Genesis[®] Malleable Penile Prosthesis is a flexible silicone elastomer device designed to be implanted into the penis for the management of erectile dysfunction.




Genesis®

The power of simplicity.

Simple implantation and use


Device construction, hydrophilic coating, and tip shape allow for a small corporotomy and better fit providing a natural aesthetic result, with unique rigidity and superior concealability.



HydroVANTAGE™ The power of choice.

With HydroVANTAGE Hydrophilic Coating

- HydroVANTAGE™ hydrophilic coating offers physician choice in selection of aqueous solution dip
- Hydrophilic coating increases the lubriciousness of the device, making it easier to implant



2-piece flexible penile prosthesis with hydrophilic coating



Trimmable

- Custom sizing to each individual’s corporal length



Easy to assemble

- No special tools required. Simply cut to size and connect the tail cap



Adjustable

- Three sizes of tail caps available for adjusting the prosthesis length after trimming

Ordering information

Genesis® Malleable Penile Prosthesis

Diameter	Length cm		Order Ref.
	mm		
		Min.	Max.
9.5		14	23
11		16	25
13		18	27



Reference
1. Habous M, Malleable (Semi-Rigid) penile prosthesis. J Sex Med. 2015 Oct;12 (10):1984-1988.