


Virtue[®]

Male Sling System

The Virtue[®] Male Sling System is a permanent, synthetic suburethral sling. It is designed for the surgical treatment of male stress urinary incontinence (SUI).



Because moderate matters.

Moderate stress incontinence deserves a precise solution

The Virtue Male Sling's dual mode of action¹ fills the gap between standard slings and invasive treatments like the AUS². With its innovative four-arm design, the Virtue sling offers anatomical support³ to improve continence outcomes.⁴

The Virtue sling difference

Prepubic arms

provide urethral **compression** for increased resistance⁶

Transobturator arms


provide proximal relocation through **elevation** of the bulbous urethra⁶

Anatomical fit³ for those in-between²

Its dual adjustment modes and four-arm design provide an anatomical fit³ that addresses the needs of moderate incontinence patients.² Our innovative quadratic design provides optimal support while reducing complications.⁴

Anatomical fit³ for ongoing sling support⁶

- 2 Dual tension adjustments¹ contributing to patient's continence
- 4 Quadratic fixation⁶ limits sling loosening and/or movement



Minimally invasive for fast recovery^{*,4}

The Virtue sling allows your patients to return to their lifestyles fast^{*,7} with its minimally invasive approach.⁴

- Strategic alternative to more invasive intervention like the AUS
- Increased quality of life⁴

Ordering Information

Product	Order number
Virtue Male Sling System 1 sling, 1 single-use introducer	50020



 Before use, always consult the 'Instruction For Use' document delivered with the products/devices.

Ostomy Care | Continence Care | Wound and Skin Care | Interventional Urology | Voice and Respiratory Care

www.coloplast.com Coloplast and the Coloplast logo are trademarks of Coloplast A/S.
© 2025 Coloplast A/S. All rights reserved. PM-37005